

# **Aerodynamic Aerosol Concentrators**

---

**Shawn Conerly  
John S. Haglund  
Andrew R. McFarland**

**Aerosol Technology Laboratory  
Texas A&M University/Texas Eng. Exp. Station**

**Presented at:  
Joint Services Scientific Conference on  
Chemical and Biological Defense Research  
Towson, MD**

**November 19, 2003**

Report Documentation Page				Form Approved OMB No. 0704-0188	
Public reporting burden for the collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to a penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.					
1. REPORT DATE <b>19 NOV 2003</b>		2. REPORT TYPE <b>N/A</b>		3. DATES COVERED <b>-</b>	
4. TITLE AND SUBTITLE <b>Aerodynamic Aerosol Concentrators</b>				5a. CONTRACT NUMBER	
				5b. GRANT NUMBER	
				5c. PROGRAM ELEMENT NUMBER	
6. AUTHOR(S)				5d. PROJECT NUMBER	
				5e. TASK NUMBER	
				5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) <b>Aerosol Technology Laboratory Texas A&amp;M University/Texas Eng. Exp. Station</b>				8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSOR/MONITOR'S ACRONYM(S)	
				11. SPONSOR/MONITOR'S REPORT NUMBER(S)	
12. DISTRIBUTION/AVAILABILITY STATEMENT <b>Approved for public release, distribution unlimited</b>					
13. SUPPLEMENTARY NOTES <b>See also ADM001851, Proceedings of the 2003 Joint Service Scientific Conference on Chemical &amp; Biological Defense Research, 17-20 November 2003. , The original document contains color images.</b>					
14. ABSTRACT					
15. SUBJECT TERMS					
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT <b>UU</b>	18. NUMBER OF PAGES <b>20</b>	19a. NAME OF RESPONSIBLE PERSON
a. REPORT <b>unclassified</b>	b. ABSTRACT <b>unclassified</b>	c. THIS PAGE <b>unclassified</b>			

# Overview

---

- **Detector Considerations**

- Response threshold on the order of  $10^4$  particles.
- Air handling system is a major power consumer in the detector system.

- **Goal:** Develop High-Throughput, Low-Power Bioaerosol Concentrators.

- **Slot-Nozzle Virtual Impactors**

- Smaller nozzle dimensions reduce power consumed for fixed cutpoint and flow rate concentrator.
- Slot nozzles are easier to manufacture for high through-put and resist fouling to a greater degree than round nozzle virtual impactors.


- **Prototype Devices**

- LSVI and CSVI Units
- Wall loss mechanisms
  - Misalignment of acceleration and receiver nozzle.
  - Receiver nozzle profile.
  - Acoustic resonance.
- Experimental Data

- **Research Plans**


AEROSOL TECHNOLOGY LABORATORY  
Texas A&M University

# Concept of Virtual Impaction



# Theoretical Power


(slot nozzle virtual impactor operating at  
500 L/min and with 0.8  $\mu\text{m}$  AD cutpoint)



Texas A&M University AEROSOL  
TECHNOLOGY LABORATORY

# Nozzle Geometry of Current Slot Nozzle Prototypes

---



Acceleration nozzle width:  
305  $\mu\text{m}$  (0.012")

Receiver nozzle width:  
457  $\mu\text{m}$  (0.018")


Slot length:  
89 mm (3.500")

Receiver nozzle taper:  
Compound, 508  $\mu\text{m}$  (0.020") step,  
60 degrees, full step

# Prototype Slot Nozzle Virtual Impactors

■ LSVI Prototype


■ CSVI Prototype



Texas A&M University AEROSOL  
TECHNOLOGY LABORATORY


# Multi-Stage LSVI Bioaerosol Concentrator

Plan and Cross-Section Views of  
2 Stage Virtual Impactor





# Multi-Stage CSVI Bioaerosol Concentrator



# Comparison of Particle Deposition Evident on Receiver Nozzles Without (a) and With (b) Secondary Expansion


---

(a)


(b)



# Numerical Prediction of Particle Tracks (1.0 micron and 5.0 micron) with 0.002" Nozzle Misalignment





&M Un  
IOLOG





# Sound Recordings of LSVI Open to Ambient Environment at Test Flow Rates

---


■ LSVI 42 L/min 

■ LSVI 47 L/min 


■ LSVI 57 L/min 

■ LSVI 66 L/min 


# Sound Power Spectrum of LSVI at 47.2 L/min




# Peak Sound Intensity of LSVI Unit



# Frequency of Peak Sound Component




# Collection Efficiency of LSVI Minor Flow




Texas A&M University AEROSOL  
TECHNOLOGY LABORATORY



# Wall Losses on Critical Zone of LSVI




# Collection Efficiency of CSVI Minor Flow




Texas A&M University AEROSOL  
TECHNOLOGY LABORATORY

# Wall Losses on Critical Zone of CSVI



Texas A&M University AEROSOL  
TECHNOLOGY LABORATORY

# LSVI and CSVI Minor Flow Collection Efficiency on Stokes Number Basis



# Research Plans

---

- **Construct Single-Stage LSVI Concentrator System**
  - Eight LSVI units operating in parallel for 300 L/min total sampling rate.
  - Design flow control instrumentation.
  - Test with monodisperse aerosols: PSL (non-viable), single-spore BG, monodisperse BG clusters.
  
- **Identify and Eliminate Acoustic Resonance in LSVI units**
  - Determine resonance source (Helmholtz, standing-wave).
  - Tapered major flow path.
  - Materials maximizing acoustic absorption.
  
- **Construct new CSVI unit**
  - 0.012" Acceleration Slot Width.
  - Incorporate resonance elimination features.
  - Add additional stages.
  - Test with monodisperse aerosols: PSL (non-viable), single-spore BG, monodisperse BG clusters.